

EN • pHormulis Labò

Device for cosmetic use. Does not perform diagnoses or treatments. CE mark

Cod. PLAE112

Distributor

Callegari Benelux BV

Description and Intended Use

Professional machine for preparing fresh, customized skincare treatments

Technical specifications

Weight	6 kg (unit only) 6.8 kg (with power supply)
Dimensions (in mm)	height: 320 mm – 326.5 mm. Width: 180 mm. Depth: 215.75 mm – 254.87 mm.
Operating conditions	Between 10 °C and 30 °C. Protect the unit from: direct UV exposure and humidity.
Storage conditions	With the water tank empty, the drip tray empty, and the purging cycle completed. At a temperature between 10 °C and 45 °C, away from dust and humidity
Power supply	100-240 V, 50-60 Hz.
Rated power	150 W (24 V x 6.25 A).
Power consumption	13.2 Wh on average per preparation.
Protection class	IP52
Connectivity	Bluetooth Low Energy v4
Display	Touchscreen 480 x 320 pixels
Water dispensing	Via peristaltic pump
Water temperature control	Heating up to 95°C.
Capsule Identification	Via QR code scanning.
Capsule Dosing	Managed by a stepper motor
Capsule Heating	Temperature control for large capsules up to 80°C.
Mixture Agitation	Agitation speed control up to 6700 rpm
Mixture cooling	Via fan
Hygiene	Automated disinfection system
Compliance	Complies with ISO 22716, verified by Bureau Veritas.
Process water	Demineralized
Capsules in the carousel	1 texture / 0, 1 o 2 active ingredients / 0 o 1 perfume